

A new approach to characterization of carbon pore structure

Jacek Jagiello* and Jeffrey Kenvin

Micromeritics Instrument Corporation, Norcross, Georgia, USA

*jacek.jagiello@micromeritics.com

For years, the characterization of carbon pore size distribution (PSD) has been dominated by the analysis of nitrogen (N_2) isotherms measured at its boiling point (77 K). From the scientific viewpoint, however, N_2 is not the most appropriate molecular probe for the PSD analysis because of its significant quadrupole moment that may influence the adsorption of N_2 molecules due to interactions with polar surface sites. For this reason, IUPAC Technical Report 2015 [1] recommended Ar as more reliable for the PSD analysis. In practice, however, using liquid argon as a cryogen for argon adsorption measurements may be challenging due to its high cost and limited availability in various parts of the world. To reconcile this problem, we have recently proposed using O_2 instead of N_2 or Ar [2]. We have chosen O_2 because its quadrupole moment is less than one-third of the value reported for N_2 . Moreover, we demonstrated a quantitative agreement between the PSD results derived from the adsorption isotherms of O_2 and N_2 measured at 77 K, and Ar at 87 K on representative carbon samples.

On the other hand the IUPAC Technical Report [1] recommended using CO_2 at 273 K for the analysis of ultramicroporous carbons, because at this temperature CO_2 diffuses faster into micropores than N_2 at 77 K. This recommendation, however, is inconsistent with the premise of using molecules with low susceptibility to specific interactions with polar sites, because the quadrupole moment of CO_2 is even larger than N_2 .

To improve the diffusion into ultramicropores, we propose in this work to use H_2 in combination with O_2 isotherms measured at 77 K. H_2 molecule has a very small quadrupole moment and smaller diameter than CO_2 . It has been shown [3] that H_2 isotherms may provide meaningful PSD results in an extended range of pore sizes when analyzed simultaneously with argon isotherms. In this work, we use 2D-NLDFT models to simultaneously analyze H_2 and O_2 isotherms measured at 77 K. We show comparisons of PSD results derived from single isotherms and the dual analysis of both isotherms. The approach introduced in this work is advantages compared to the earlier presented “dual analysis method”[4] based on N_2 and CO_2 isotherms, because of the lower quadrupole moments of O_2 and H_2 . It is also more practically convenient because of the same cryogen, liquid N_2 , used for both analyses.

[1] M. Thommes, K. Kaneko, A.V. Neimark, J.P. Olivier, F. Rodriguez-Reinoso, J. Rouquerol, K.S. Sing, Physisorption of gases, with special reference to the evaluation of surface area and pore size distribution (IUPAC Technical Report), *Pure and Applied Chemistry* 87 (2015) 1051-1069.

[2] J. Jagiello, J. Kenvin, Consistency of carbon nanopore characteristics derived from adsorption of simple gases and 2D-NLDFT models. Advantages of using adsorption isotherms of oxygen (O_2) at 77 K, *Journal of Colloid and Interface Science* 542 (2019) 151-158.

[3] J. Jagiello, W. Betz, Characterization of pore structure of carbon molecular sieves using DFT analysis of Ar and H_2 adsorption data, *Microporous and Mesoporous Materials* 108 (2008) 117-122.

[4] J. Jagiello, C. Ania, J.B. Parra, C. Cook, Dual gas analysis of microporous carbons using 2D-NLDFT heterogeneous surface model and combined adsorption data of N_2 and CO_2 , *Carbon* 91 (2015) 330-337.